



Administration Schedule 2

My Container is:

Order code:

- ☐ Flocare Container 0.5L 40441
- ☐ Flocare Container 1L 40455
- ☐ Flocare Top Fill Reservoir 1.3L 70139

Accessories:

Order code:

- ☐ Flocare Two Pack Connector 69915
- ☐ Flocare ENFit Bolus Adaptor 89740

Syringes:

Feeding Tube:

You will have one of the following feeding tubes:

- 1. Nasogastric** – often referred to as an NG tube. This is a narrow plastic tube, which is inserted through your nose and passes down into your stomach
- 2. Percutaneous Endoscopic Gastrostomy** – often referred to as a PEG. This is a silicone / polyurethane tube which is inserted directly into your stomach through the skin
- 3. Jejunostomy** – often referred to as JEJ. This is a silicone / polyurethane tube which is inserted directly into your small intestine through the skin
- 4. Button** – A low profile feeding tube that must be used with an extension set
- 5. Naso-Jejunostomy** – A feeding tube that is inserted through the nose into the intestine
- 6. Radiologically Inserted Gastrostomy** – Often referred to as a RIG tube. Using x-ray and scanning equipment, this tube is inserted through the skin and into the stomach

My tube is a NG / PEG / JEJ / Button / NJ / RIG tube (circle as appropriate):

Brand name / type

Size

Y connector (if applicable)

Code

Button extension set (if applicable)

Code